

# PALL Kleenpak™ Medical Filters



## Disposable Tap & Showerhead Filters

Where a positive Legionella result has been detected, rather than having to shut down the system, installing Pall medical filters allow the system to stay online and user safety is uncompromised. Pall medical filters allow vital potable water systems to remain online whilst remedial actions are planned to combat the legionella positive.

### Features:

- Immediate removal of Legionella spp for up to 1 month
- Unique filtration media engineering
- Bacteriostatic housing to minimise the risk of retrograde contamination

### Benefits:

- Instant & Cost effective protection against Legionella spp.
- Easy to install
- Long filter life, high volume throughput
- Compatible with and complimentary to thermal & chemical systemic treatments

## Product Specification:

Filter Life	Maximum Operating Pressure	Maximum Continuous Influent Temperature	Maximum Temperature Exposure	Membrane Rating
31 Days	5 bar	60°C	70°C for a total cumulative period of 30 minutes over the life of the filter	0.2 µm Supor, incorporating pre-filtration layer

## Connections & Adaptors (Various available to suit all connection types)



EPSCO Ltd: Tel: +44(0)1738 446363 Fax: +44(0)1738 442265 E-mail: [info@epsco.co.uk](mailto:info@epsco.co.uk)!