

3M Clean-Trace™ NG Luminometer ATP Tester



3M Clean-Trace™ NG Luminometer ATP Tester

Hand held, lightweight NG Luminometer to measure ATP tests and to capture data for later analysis. ATP surface and water tests are used to determine the level of contamination in a sample, results are available in seconds.

Ideal for monitoring hygiene and risks in critical food and production systems.

- 1. Obtain The Sample** - Accurate, easy to use hand held testing with results in just 30 seconds.
- 2. Get Reliable Results** - Highest degree of repeatability, giving you the lowest variability and most accurate data on the market.
- 3. Track Your Trends** - Data trending software can be optimised to identify workplace trends.



You may also need...



3M™ Clean-Trace™ Surface ATP

Unique liquid stable enzyme, the test offers superior levels of repeatability. Real-time ATP hygiene test which gives an objective indication of cleanliness in seconds (Pack of 100).

EPSCO Ltd: Tel: +44(0)1738 446363 Fax: +44(0)1738 442265 E-mail: info@epsco.co.uk